





# High performance sterilisation paper Steripap

In the healthcare every precaution is very important to ensure risk prevention. For this reason we work rigorously to offer the most effective and safely products that provide the quality expected by our customers.

Dorsan® Steripap is a high purity paper made with alphacellulose fibers, free from all environmental contaminants. Used in medical field for wrapping instruments which are to be sterilized by steam (121-134°C) in Autoclaves, by means of ethylene oxide, by Gamma irradiation (50 kGy) or formaldehyde.

We offer a wide range of sizes in colors white, blue or green. The quality of our product is subject to well-defined and regularly monitored standards by means of the standards DIN EN ISO 9001.

### **Features**

Easy wrapping and packaging of instruments and medical devices.

They allow the contents of the packages remain totally sterile.

Protect the contents of the package during storage and transport and promote the aseptic presentation at point of use.

## **Available Colours**

Green

Blue

White









### **Sterilisation paper Steripap specifications**

| Size         | Quality                                                  | Packaging options                                                                                                |
|--------------|----------------------------------------------------------|------------------------------------------------------------------------------------------------------------------|
| 40 x 40 cm   | 58 g/m² blue, green, white added sizing agent, wood-free | 100 to 500 sheets per<br>pack according to<br>sheet size, shrink-<br>wrapped in a tray or in<br>an outer carton. |
| 50 x 50 cm   |                                                          |                                                                                                                  |
| 60 x 60 cm   |                                                          |                                                                                                                  |
| 75 x 75 cm   |                                                          |                                                                                                                  |
| 90 x 90 cm   |                                                          |                                                                                                                  |
| 100 x 100 cm |                                                          |                                                                                                                  |
| 120 x 120 cm |                                                          |                                                                                                                  |

Please ask for a bespoke product offer.

# Recommendations

Store the product in a dry place, free of dust and with little sunlight at a temperature of 10 to 20 ° C.

Use sequential wrapping system to obtain the best barrier properties and more aseptic presentation.























All indicated data to be understood as typical average values. This product is manufactured in conformance with the Quality Management standard DIN EN ISO 9001. We reserve the right to change these informations without any previous notice. v01 © Dorsan Filtration

